Mepilex® Border Flex

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**Références 1.** Nelson, D. A New Bordered Foam Dressing Technology Improves Wound Outcomes and Satisfaction While Reducing Dressing Utilization in Acute Care. **Poster présenté au congrès annuel de la WOCN**, Nashville, 2019 2. Tyson, L. Study First: Driving the Case for Improving Hospital Wound Care. **Poster présenté au congrès annuel de la WOCN**, Nashville, 2019 3. Mölnlycke Health Care. Mepilex® Border Flex – Estimate of spreading area using dot pattern on backing film. Rapport n° PD-528872. Données internes. 4. Mölnlycke Health Care. Mepilex® Border Flex Product Manual – Conformability PD-528870. Données internes. 5. Mölnlycke Health Care. Mepilex® Border Flex-Bacteria encapsulation. Rapport n° PD-537072. Données internes. 6. Mölnlycke Health Care. Mepilex® Border Flex Product Manual – Spreading layer PD-528874. Données internes. 7. Mölnlycke Health Care. Product Manual – Maintained dressing profile. Rapport n° 525458. Données internes. 8. White, R. A multinational survey of the assessment of pain when removing dressings. Wounds UK 2008;4(1):14-22. 9. Mölnlycke Health Care. Mepilex® Border Flex – Absorption and dispersion tests on inclined plane, viscous test solution. Rapport n° PD-528871. Données internes. 10. Mölnlycke Health Care. Product manual Initial abs different test solutions. Rapport n° PD-535352. Données internes. 11. Mölnlycke Health Care. Mepilex® Border Flex – Fluid handling capacity compared to Urgotul Foam Absorb Border, Aquacel Foam, Allevyn Classic Gentle Border, Allevyn Life, Optifoam Gentle, KerroFoam Gentle Border. Rapport n° PD-527642. Données internes. 12. Mölnlycke Health Care. Product Manual Mepilex® Border Flex – Viral penetration. Rapport n° PD-535090. Données internes. 13. Opinion Matters 'Professional Health Care Treatment for Wounds' survey conducted June-July 2018 in UK. **Données internes 14.** Woo, K. & Bergström, C. An Open Randomised Cross-over Investigation Assessing Perceived Pain at Dressing Change Comparing A Soft Silicone Dressing, Mepilex® Border, with a Adhesive hydrocellular polyurethane dressing, Allevyn® Adhesive in Patients with chronic ulcer (Clinical Investigation Report (MXB408, Part A, Mölnlycke Health Care), 2007). 15. Van De Looerbosch, D., Anderson, A., Husmark, A. & Körner, V. An open randomised explorative study comparing Mepilex® Border with Telle in patients with pressure ulcers according to stage II. EPUAP guidelines (Clinical Investigation Report (MPX/B201), Mölnlycke Health Care, 2001). 16. Mölnlycke Health Care. Mepilex® Border Flex – Waterproofness. Rapport n° PD-532095. Données internes.

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Les marques, noms et logos Mölnlycke, Mepilex et Safetac sont déposés par une ou plusieurs sociétés du groupe Mölnlycke Health Care pour le monde entier. ©2020 Mölnlycke Health Care AB. Tous droits réservés. BF-WC00032002